

EPSA[®] *safety first*

**ARMFLEX LARYNGECTOMY
TUBE WITH CUFF**

ARMFLEX LARYNGECTOMY TUBE WITH CUFF

CODE: 9220

1 PVC tube reinforced with flexible stainless-steel spiral, polished distal orifice, atraumatic. Low-pressure cuff, pilot balloon and retention valve. Standard connector.



Available Sizes

4,0 / 5,0 / 6,0 / 7,0 / 8,0 / 9,0 / 10,0 mm

STERILE, SINGLE PACK, SINGLE USE.

Recommendation for use

Recommendation of a professional on the use of the technology.

The Armflex Laryngectomy tube with cuff, is indicated in patients undergoing laryngeal surgery.

The cuff seals the tracheal lumen, allowing adequate control of gas movement during mechanical ventilation and preventing the passage of fluids and secretions from the upper airway into the bronchial tree.

The wall of the cannula is reinforced with a stainless-steel spiral that provides high resistance to kinking and the collapse of its lumen by extrinsic compression.

