

**EPSA<sup>®</sup>** *safety first*

ESCHMANN INTUBATION STYLETS

# ESCHMANN INTUBATION STYLETS

CODE: 9230


1 PVC stilet graduated in cm. with curved end.

**STERILE, SINGLE PACK, SINGLE USE.**

## Recommendation for use

Recommendation of a professional on the use of the technology.

The Eschmann Intubation Stylet is a special auxiliary device to facilitate the procedure in case of difficult endotracheal intubation. It has a slightly curved distal end, which allows, in the blind introduction, to feel the contact of the tip of the same with the tracheal rings and confirm its position within the airway.



### Available Sizes

2,6 mm x 700 mm.  
5,2 mm. x 700 mm.