



EPSA[®] *safety first*

**TRACHEOSTOMY CANNULA
WITHOUT CUFF**

TRACHEOSTOMY CANNULA WITHOUT CUFF

CODE: 9751

1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic. Flexible fixation flange with standard connector.

1 insertion aid with olive tip.
1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.



Available Sizes

5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 /
8,5 / 9,0 / 9,5 / 10,0 mm

Recommendation for use

Recommendation of a professional on the use of the technology.

The Tracheostomy Cannula without cuff is indicated in tracheostomized patients to allow the maintenance of respiratory function, with adequate management of the airway, in relation to hydration, aspiration of secretions and endoscopic maneuvers for diagnosis or treatment.

